

NEW!

A-dec 200 with
Traditional Delivery,
Chair Mount Light,
Assistant's Instrumentation,
Cuspidor & Fibre Optics

\$20,155*

A-dec 200

Designed for success

OPTIONS: Intra-Oral Camera \$1,000
Air Motor \$1,000
Scaler \$1,400



* Additional charge
for installation.

a dec[®]
reliablecreativesolutions

SALES - SERVICE - SUPPORT

BRISBANE HEAD OFFICE

Showroom P 07 3862 6288
Unit 6, 43-49 Sandgate Rd F 07 3862 3109
ALBION QLD 4010 E sales@dentaldepot.com.au

MELBOURNE OFFICE


Showroom P 03 9676 2700
202 Turner Street F 03 9676 2005
PORT MELBOURNE VIC 3207 E sales@dentaldepot.com.au

ALL MAIL

Postal
PO Box 788,
ALBION QLD 4010



FREE iPod Touch with every Sonoswiss Ultrasonic Bath



Useful and practical accessories. The cover can also serve as a drip-off collecting basin for the wet basket.



Connection for drain duct at the back (not on 3Ltr model). The workplace can be chosen with or without a wash basin.



The housing, made completely of Stainless Steel, allows hygienic cleaning and disinfection of the unit surfaces.



Protecting rim to keep the operating panel free of cleaning medium. Any spilled water is guided past the operating panel.

The new Sonoswiss removes organic contaminants, sterilizing residues, silicate stains and other substances. Its sweep function ensures that instruments are cleaned reliably no matter where they are placed in the bath. With built in heater, grease and wax are removed when heated above their fusion temperature.

Available in...

SW 6H	(6 Litre)	\$1,600
SW 12H	(12 Litre)	\$2,500



Plant Room Specials... COMPRESSORS



PERFORMANCE and RELIABILITY

HIGH PERFORMANCE COMPRESSORS COVERED BY A 5YR WARRANTY.

It's widely understood that dental clinics need an oil free compressor. However, what isn't so well understood is that not all oil free compressors are the same.

WHAT ARE THE DIFFERENCES?

Two simple things: performance and reliability. Concerning performance, we recommend you take the time to compare levels of air delivery between manufacturers. And as for reliability, well, we think the 5 year warranty says a lot.

If you'd like to know more about dental compressed air systems and their effect on your clinic, talk to DENTAL DEPOT... because we're the specialists.

CAT-K100E	(1 Surgery)	\$5,500
CAT-K200E	(2-3 Surgery)	\$7,300
CAT-K300E	(3-5 Surgery)	\$8,300



TURBOSMART – Intelligent Suction

SMARTER THAN FIXED-SPEED...

Turbo SMART isn't like other suction systems: it varies its speed according to your surgery requirements. Turbo SMART's motor adapts to provide you with consistent high airflow and vacuum, regardless of whether other surgeries are using the suction unit.

No more loss of power!

ROOM TO GROW...

Two Turbo Smart models are available, one catering for two surgeries operating simul-taneously, and the other up to four surgeries operating simultaneously. Both models are, in fact, identical machines containing different software.

This means that as your number of surgeries expands, you can update your model by simply entering a new password!

CAT-C1A	(2-3 Surgery)	\$7,800
CAT-C2A	(3-4 Surgery)	\$9,500

Even less without Amalgam Separation System! 2 YEAR WARRANTY



ALL PRICES
QUOTED
ARE GST
INCLUSIVE.

Safety,
purely and simply



MAY / JUNE SPECIALS

Buy the
Quaz Autoclave
and Traceability
System for
\$11,750



EXCLUSIVE TO DENTAL DEPOT*

2 years or 2000 cycles
warranty.

European made.

Serviced via the
Dental Depot Group

* Qld, Vic, Tas, SA, & NT



Label Data includes:

- Doctor's Name
- Quaz serial number
- Date and time of the cycle
- Cycle number
- Program type
- Barcode
- Expiry date.



ALL PRICES
QUOTED
ARE GST
INCLUSIVE.

Ti-Max X Series START-UP SPECIAL

HANDPIECE PACKAGE Solid Titanium

- Solid Titanium Body with Scratch Resistant DURACOAT
- Cellular Glass Optics • Clean Head System • Push Button Chuck

PACKAGE INCLUDES:

3 x Ti-Max X High Speed Optic Handpieces

- Choice of any three: X500L mini head
X600L standard head
X700L torque head

1 x NSK LED Coupling

3 x Ti-Max X Contra-angle Handpieces

- Choice of any three: X95L 1:5 speed increasing handpiece
X25L 1:1 direct drive handpiece
X15L 4:1 speed reducing handpiece
X65L 1:1 straight handpiece
X57L 16:1 speed reducing prophylaxis handpiece

2 YEAR WARRANTY 



OUR PRICE
\$8,500
WAS \$9,999

S-Max M Series START-UP SPECIAL

HANDPIECE PACKAGE Stainless Steel

- Stainless Steel Body • Cellular Glass Optics
- Clean Head System • Push Button Chuck

PACKAGE INCLUDES:

3 x S-Max High Speed Optic Handpieces

- Choice of any three: M500L mini head
M600L standard head

1 x NSK LED Coupling

3 x S-Max X Low Speed Handpieces

- Choice of any three: M95L 1:5 speed increasing handpiece
M25L 1:1 direct drive handpiece
M15L 4:1 speed reducing handpiece
M65 1:1 straight handpiece

1 YEAR WARRANTY 



OUR PRICE
\$5,695
WAS \$6,699

ADD-ON PACKAGE 1:

CLINICAL MICROMOTOR and ULTRASONIC SCALER

PACKAGE INCLUDES:

1 x NLX Plus integrated electric e-type micromotor with Multipad

- Easily integrated into your existing chair
- Endodontic function
- Torque reaching 4.0 Ncm
- Wide rotational range of 100-40,000 min-1
- Auto-reverse on/off and stop functions

1 x NSK Varios 170 integrated multifunction ultrasonic scaler with LED illumination

- Choice of over 70 tips
- Lightweight handpiece with twin LED lights
- Powered by NSK's new iPiezo engine



OUR PRICE
\$4,760
WAS \$5,599

Multipad touch panel controls the NLX Pls and Varios 170 for ease of use and complete user control

ADD-ON PACKAGE 2:

OPTIC AIR MOTOR, AIR POLISHER and POWDER

PACKAGE INCLUDES:

1 x Ti205L optic air motor

- Speed: max 20,000 RPM
- Internal water spray

1 x Prophy-Mate Neo air polisher

- Ultra lightweight and portable (transport case included)
- 360° easy swivel for stress free and comfortable use
- Connectivity with all major coupling systems
- 60° and 80° nozzles included for easy access

4 x bottles of Flash Pearl air polishing powder (300g ea)

- Calcium based spherical particles



OUR PRICE
\$2,805
WAS \$3,299

ADD-ON PACKAGE 3:

HANDPIECE MAINTENANCE

CHOICE OF ONE:

1 x Care3 Plus C1 PLUS 2 x PTL adaptors

One rotation e-type connection and two high speed connections

1 x Care3 Plus C2 PLUS 1 x PTL adaptor

Two rotation e-type connections and one high speed connection

- A fail safe maintenance process
- Easily cleaned
- Very user friendly – simple to operate
- Rotates handpieces during the cycle
- Economical, efficient and time saving, 2 minute cycle



OUR PRICE
\$3,395
WAS \$3,999

CARE 3 PLUS C2



ALL PRICES QUOTED ARE GST INCLUSIVE.

THERMAL DISINFECTOR – Protection, Economy & Validation

Processing instruments in the Miele washer-disinfector cleans and disinfests instruments in a closed system. The process is easy to validate and revalidate.

With its fungicidal, bactericidal and virus-deactivating properties, thermal disinfection achieves the highest possible level of protection against risk of infection.

It is only by using this process that the most essential criteria for professional instrument processing in dental practices are fulfilled. The Miele G 7881 Dental Washer-Disinfector is the result of extensive R&D in close cooperation with the dental industry.

It has been developed specifically for the cleaning of dental instruments and accessories.



WHY THE MIELE G 7881?


1. Delivers a high level clean: Exceeds OSHA standards for Microbial contamination reduction of 99.9%.
2. No chemical bath: Disinfection Vario Program delivers thermal disinfection in the final rinse when temperature reaches 90° C and holds for 10 mins.
3. The Right Capacity: A dual level chamber can accommodate baskets, inserts and mesh trays. Cleans hinged instruments and cassettes.
4. Simplifies reprocessing: Offers 4 steps in 1 - Pre-soak, Rinse, Clean and Dry.
5. Interchangeable baskets available: For customization including General Dentistry and Oral Surgery packages.
6. Simple process: Load machine, select a cycle, and let the machine do the rest. Eliminates pre-soaking, hand scrubbing, rinsing and drying.
7. Prepares instruments for sterilization more quickly and safely.



G 7881 AE
Thermal Disinfector with stainless steel exterior.



U 874
Standard lower basket. Empty, for various inserts.




O 190
Standard upper basket. Empty, for various inserts.



E 146 Mesh tray (1/6) with lid for hand instruments. 25/32" x 57/8" x 827/32", constructed in 1/16" and 3/32" wire mesh.



E 363 Mesh tray (1/6) As E 146, but constructed of 1/32" wire mesh.



E 378 Full universal mesh tray (1/1) 315/16" x 183/32" x 183/32".



E 379 Half universal mesh tray (1/2) 315/16" x 71/16" x 183/32" for various utensils, made of 1/16" wire mesh.



E 328 Support insert Suitable for use in Mesh trays E146 and E363. For upright support of single and double ended instruments.



E 417 Insert (2/5) for lower basket, for miscellaneous items.



E 131/1 Insert for lower basket only, suitable for 5 cassettes.



E 130 Insert (1/2) for lower basket, for 10 trays.



E 339 Insert (3/5) for lower basket, for 16 tray bases.



E 373 Mesh tray (1/6) as E 363, but with instrument holders.



E 338 Insert (3/5) for upper and lower baskets, for 8 half trays.



E 523 Insert Holds 6+ cassettes. Fits in upper or lower-basket. Recommended insert for Hu-Friedy cassettes.

Miele is known throughout the world for high quality, high performance aqueous cleaning systems.

For more than thirty years, Miele has provided dentists, hospitals and laboratories with superior cleaning systems.

Our resulting expertise, combined with a commitment to unsurpassed quality and performance, is evident in the basic features found in every Miele Cleaning System.

* Available in Queensland only

Miele
PROFESSIONAL

Are your X-rays due for compliance testing?

Dental Depot is fully accredited for
X-ray Compliance Testing.

We use state-of-the-art radiation
detection equipment providing accurate
testing with a minimum of disruption to
your practice whilst adhering strictly to the
relevant State and National regulations
and guidelines.

Customers receive an Assessment Report
including a Certificate of Compliance.

Dental Depot offers this service at
\$275 per intra-oral x-ray.

If you decide to purchase an
intra-oral x-ray from us, we will include
the compliance test for \$110!

Once you become a valued client of
Dental Depot for x-ray compliance testing,
you will automatically be sent a reminder
when the test becomes due again.

So, if you are in doubt about whether your
X-ray equipment complies,
or you would like further information
on the services we provide, give our
friendly staff a call.

We Service What We Sell

DIGITAL X-RAY PACKAGE Includes KODAK 2100 X-Ray, Size 1 Digital Sensor, Installation and Compliance Test



\$13,945

The Kodak 2100 Intraoral
X-Ray System is specifically
designed for digital systems.
It offers high frequency
technology, ease of use and
reliability. It combines high quality
results for accurate diagnosis with a flexible
and smart design.

Featuring a high frequency generator to help obtain
high image quality while minimizing radiation exposure,
by up to 30 percent, the Kodak 2100 system is compatible
to every type of imaging technology.

The system operates with both film and digital
radiography, including the Kodak CR 7400 and RVG digital
radiography systems. The 2100 system includes a timer
which can be mounted remotely or directly on the unit.

Kodak Dental Systems

MAY / JUNE SPECIALS

ALL PRICES
QUOTED
ARE GST
INCLUSIVE.

NEW MODEL!
DD-4000 SERIES SURGICAL MICROSCOPE
with cool, white LED illumination



Standard Beam Splitter and SLR
Camera Adaptor



3D Rotation Binocular Assistant Scope



Easily maneuverable
optical head with
tension control knobs
and two handles



**Free iPad
with every
Surgical
Microscope
purchase!**

LED Version
\$13,500

The DD-4000 series Surgical Microscope is based on the DD-3000 Series, but comes with 6-step magnification and LED illumination. It offers outstanding optical performance with brilliant apo-achromatic optics and excellent ergonomics. It combines high performance and reliability.

The ergonomic design of the microscope gives precise, fatigue-free operation with brilliant 3 dimensional images at a large depth of focus.

It is suitable to be used as a Dental and ENT surgical microscope especially.

Features:

- 20 Watt LED Illumination
50000 Min. Lux
- Mobile stand
- Over 50,000 hours of LED life

Options

- Floor, wall and ceiling mounts avail.
- 3D Rotation Binocular Assistant Scope available.
- SLR digital Camera adapter



ALL PRICES
QUOTED
ARE GST
INCLUSIVE.

SALES - SERVICE - SUPPORT

BRISBANE HEAD OFFICE

Showroom P 07 3862 6288
Unit 6, 43-49 Sandgate Rd F 07 3862 3109
ALBION QLD 4010 E sales@dentaldepot.com.au

MELBOURNE OFFICE

Showroom P 03 9676 2700
202 Turner Street F 03 9676 2005
PORT MELBOURNE VIC 3207 E sales@dentaldepot.com.au

ALL MAIL

Postal
PO Box 788,
ALBION QLD 4010